

HLA-DRB5 Polyclonal Antibody



Catalog Number: 161232

Product name

HLA-DRB5 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to HLA-DRB5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 30-227 of human HLA-DRB5 (NP_002116.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:1000 - 1:2000

Validations



Western blot - HLA-DRB5 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HLA-DRB5 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.