

CCL26 Polyclonal Antibody



Catalog Number: 165232

Product name

CCL26 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CCL26

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-94 of human CCL26 (NP_006063.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - CCL26 Polyclonal Antibody


Western blot analysis of extracts of mouse testis, using CCL26 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - CCL26 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using CCL26 antibody at dilution of 1:200



(40x lens).


CCL26 Polyclonal Antibody



Catalog Number: 165232

Immunohistochemistry - CCL26 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using CCL26 antibody at dilution



of 1:200 (40x lens).

Immunohistochemistry - CCL26 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using CCL26 antibody at dilution of 1:200 (40x lens).