

# Anti-FACL4 Rabbit antibody



Catalog Number: 167232

## Product name

Anti-FACL4 Rabbit antibody

IHC: 1/50

ICC/IF: 1/50

## Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20


## Antibody description

FACL4 Rabbit polyclonal antibody

## Validations

## Preparation

Antigen: A synthetic peptide of human FACL4



## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of FACL4 in K562, C6, 3T3, HeLa cell lysates using FACL4 Rabbit pAb(1:1000 diluted). Predicted band size: 79kDa. Observed band size: 79kDa.

## Clonality


Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP



## Dilutions

WB: 1/5000

Immunohistochemistry of FACL4 in paraffin-embedded Human lung cancer tissue using FACL4 Rabbit pAb at dilution 1/50