

Anti-HCE Rabbit antibody



Catalog Number: 167732

Product name

Anti-HCE Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

HCE Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human HCE

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/100

IP: 1/20

Validations



Western blot detection of HCE in HeLa, A549, HL-60, U2OS, C6 cell lysates using HCE Rabbit mAb (1:1000 diluted). Predicted band size: 69 kDa. Observed band size: 69 kDa.