

LAP3 antibody



Catalog Number: 112142

Product name

LAP3 antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LAP3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, mouse heart tissue, SH-SY5Y cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 56kd

Preparation

This antibody was obtained by immunization of LAP3 recombinant protein (Accession Number: NM_015907). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112142(LAP3 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112142(LAP3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112142(LAP3 antibody) at dilution of 1:50 (under 40x lens)