

SFI1 antibody



Catalog Number: 115142

Product name

SFI1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SFI1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells, Hela cells. Positive WB detected in A549 cells, HepG2 cells, mouse lung tissue. Observed molecular weight by Western-blot: 148kd

Preparation

This antibody was obtained by immunization of SFI1 recombinant protein (Accession Number: BC033613). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:115142(SFI1 antibody) at dilution of 1:400



Immunofluorescent analysis of HepG2 cells using Catalog No:115142(SFI1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)