

TIMP-4 antibody



Catalog Number: 116142

Product name

TIMP-4 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

IF: 1:10-1:100

Antibody description

TIMP-4 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, BxPC-3 cells, HT-1080 cells, human heart tissue. Positive IHC detected in human gliomas tissue, human heart tissue. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of TIMP-4 recombinant protein (Accession Number: NM_003256). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IHC


Dilutions

Recommended Dilution:


Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116142(TIMP4 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:116142(TIMP4 Antibody) at dilution of 1:50 (under 10x lens)




TIMP-4 antibody



Catalog Number: 116142

Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:116142(TIMP4 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of MCF-7 cells using Catalog No:116142(TIMP4 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG