

# Ozuriftamab reference antibody



Catalog Number: 177142

## Product name

Ozuriftamab reference antibody

## Specificity

Human ROR2

## Antibody description

Anti-ROR2 Antibody (Ozuriftamab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-ROR2 Antibody (Ozuriftamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-ROR2 Antibody (Ozuriftamab) is 99.73%, determined by SEC-HPLC.



## Bioactivity: ELISA

Immobilized human ROR2 His at 2 µg/mL can bind Anti-ROR2 Antibody (Ozuriftamab), EC50=0.114 µg/mL.