

# Recombinant JAM3 protein (His tag)



Catalog Number: 548142

---

## Gene Name Synonym

date of receipt at -70°C

JAM-C; JAM-2; Junctional adhesion molecule 3;  
JAM-3; sJAM-C

## Molecular Mass

Predictes a molecular mass of 25.58 kDa. In SDS-PAGE under reducing conditions, it migrates as an approximately 27 kDa band.

## Protein Construction

DNA encoding 32-241 aa (Q9BX67) were fused with 6His tag.

## Formulation

Lyophilized from sterile PBS, pH 7.4

## Organism

Human

## Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

## Expression Host

E. coli

## Activity

Not tested.

## Reconstitution

Reconstitute at 0.25 µg/µl in sterile water for short-term storage. Reconstitution with 50% glycerol solution is recommended for longer term storage (see Stability and Storage for more details). If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used). Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

## Purity

>85%, by SDS-PAGE with Coomassie Brilliant Blue staining.

## Endotoxin

Please contact the lab for more information.

## Stability

Samples are stable for up to twelve months from