

Anti-NeuN Rabbit antibody



Catalog Number: 166942

Product name

WB: 1/1000

Anti-NeuN Rabbit antibody

ICC/IF: 1/100

Specificity

FC: 1/20

Human, Mouse, Rat

Antibody description

NeuN Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human NeuN

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, ICC/IF, FC

Dilutions

Validations



Western blot detection of NeuN in Mouse brain lysates using NeuN Rabbit pAb(1:1000 diluted). Predicted band size: 34 kDa. Observed



band size: 46-55 kDa.

Western blot detection of NeuN in Rat Brain lysates using NeuN Rabbit pAb(1:1000 diluted). Predicted band size: 34 kDa. Observed band size: 46-55 kDa.