

C11orf67 antibody

Catalog Number: 108642

Product name

C11orf67 antibody

Applications

ELISA, IHC

Immunogen

[Human C11orf67 Recombinant protein \(GST tag & His tag\)](#)

Dilutions

Recommended Dilution:
IHC: 1:20-1:200

Specificity

human, mouse, rat; other species not tested.

Validations

Antibody description

C11orf67 Rabbit Polyclonal antibody. Positive IHC detected in human ovary.

Preparation

This antibody was obtained by immunization of C11orf67 recombinant protein (Accession Number: NM_001392038). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage


Store at -20°C. DO NOT ALIQUOT

Clonality


Polyclonal

Ig Type

Rabbit IgG



Immunohistochemical of paraffin-embedded human ovary using Catalog No:108642(C11orf67 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human ovary using Catalog No:108642(C11orf67 antibody) at dilution of 1:50 (under 40x lens)