

Anti-TPPP Rabbit antibody



Catalog Number: 167052

Product name

Anti-TPPP Rabbit antibody

WB: 1/2000

IHC: 1/20

Specificity

Human, Mouse, Rat

ICC/IF: 1/20

FC: 1/20

Antibody description

TPPP Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Tubulin Polymerization Promoting Protein

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC

Dilutions

Validations



Immunohistochemistry of Tubulin Polymerization Promoting Protein in paraffin-embedded Human lung cancer tissue using Tubulin Polymerization Promoting Protein Rabbit pAb at dilution 1/20



Western blot detection of TPPP in Rat Brain lysates using TPPP Rabbit pAb(1:1000 diluted). Predicted band size: 24kDa. Observed band size: 24kDa.