

# Anti-Rsk 2/MAPKAP Kinase 1b Rabbit antibody



Catalog Number: 167852

## Product name

Rabbit IgG

Anti-Rsk 2/MAPKAP Kinase 1b Rabbit antibody

## Applications

## Specificity

WB, IHC-P, FC, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/2000

Rsk 2/MAPKAP Kinase 1b Rabbit monoclonal antibody

IHC: 1/20

FC: 1/20

## Preparation

IP: 1/20

Antigen: A synthetic peptide of human Rsk 2/MAPKAP Kinase 1b

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Western blot detection of Rsk 2/MAPKAP Kinase 1b in CHO-K1,Ramos using Rsk 2/MAPKAP Kinase 1b Rabbit mAb(1:1000 diluted)

