

TFF1 Polyclonal Antibody



Catalog Number: 163552

Product name

TFF1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TFF1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-84 of human TFF1 (NP_003216.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000


IHC 1:50 - 1:200

Validations



Western blot - TFF1 Polyclonal Antibody


Western blot analysis of extracts of BT-474 cells, using TFF1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - TFF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human mammary cancer using TFF1 antibody at dilution of 1:200 (40x lens).




TFF1 Polyclonal Antibody



Catalog Number: 163552


Immunohistochemistry - TFF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using TFF1 antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - TFF1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse brain using TFF1 antibody at dilution of



1:200 (40x lens).

Immunohistochemistry - TFF1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat lung using TFF1 antibody at dilution of 1:200 (40x



lens).


Immunohistochemistry - TFF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using TFF1 antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - TFF1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat kidney using TFF1 antibody at dilution of 1:200



(40x lens).

Immunohistochemistry - TFF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using TFF1 antibody at dilution of 1:200 (40x lens).



TFF1 Polyclonal Antibody



Catalog Number: 163552

Immunohistochemistry - TFF1 Polyclonal Antibody human stomach using TFF1 antibody at dilution of 1:200 (40x lens).
Immunohistochemistry of paraffin-embedded