

# ZNF169 antibody



Catalog Number: 116952

## Product name

ZNF169 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF169 Rabbit Polyclonal antibody. Positive WB detected in mouse small intestine tissue, HeLa cells, HepG2 cells, mouse liver tissue. Observed molecular weight by Western-blot: 68 kDa

## Preparation

This antibody was obtained by immunization of ZNF169 recombinant protein (Accession Number: XM\_017014364). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:116952(ZNF169 Antibody) at dilution of 1:600