

# EIF3M antibody



Catalog Number: 107252

## Product name

EIF3M antibody

## Immunogen

[Human EIF3M Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

EIF3M Mouse Monoclonal antibody. Positive WB detected in Hela cells, HepG2 cells, MCF7 cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 43kD

## Preparation

This antibody was obtained by immunization of EIF3M recombinant protein (Accession Number: NM\_006360). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107252(EIF3M antibody) at dilution of 1:1000



Immunofluorescent analysis of MCF-7 cells, using EIF3M antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red).