

ARFGAP3 antibody



Catalog Number: 108252

Product name

ARFGAP3 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARFGAP3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive IF detected in A431 cells, HepG2 cells. Positive WB detected in HepG2 cells, HeLa cells, Jurkat cells, mouse brain tissue. Observed molecular weight by Western-blot: 52-55kd, 57-66kd

Preparation

This antibody was obtained by immunization of ARFGAP3 recombinant protein (Accession Number: NM_014570). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108252(ARFGAP3 antibody) at dilution of 1:2000




Immunofluorescent analysis of A431 cells using Catalog No:108252(ARFGAP3 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



ARFGAP3 antibody

Catalog Number: 108252

Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:108252(ARFGAP3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:108252(ARFGAP3 Antibody) at dilution of 1:50 (under 40x lens)