

# SPTBN1 antibody



Catalog Number: 115652

## Product name

SPTBN1 antibody

## Immunogen

[Human SPTBN1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SPTBN1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse heart tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 275 kDa

## Preparation

This antibody was obtained by immunization of SPTBN1 recombinant protein (Accession Number: NM\_178313). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115652(SPTBN1 Antibody) at dilution of 1:600



Immunofluorescent analysis of HEK-293 cells using Catalog No:115652(SPTBN1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)