

LSM8 antibody



Catalog Number: 112352

Product name

LSM8 antibody

Immunogen

[Human LSM8 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LSM8 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HL-60 cells. Observed molecular weight by Western-blot: 11 kDa

Preparation

This antibody was obtained by immunization of LSM8 recombinant protein (Accession Number: NM_016200). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112352(LSM8 Antibody) at dilution of 1:300