

Anti-NUDT9 Rabbit antibody



Catalog Number: 167062

Product name

Anti-NUDT9 Rabbit antibody

Specificity

Human, Rat

Antibody description

NUDT9 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human NUDT9

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/5000

IHC: 1/50


ICC/IF: 1/50

FC: 1/20

Validations



Western blot detection of NUDT9 in K562, HeLa cell lysates using NUDT9 Rabbit pAb(1:1000 diluted). Predicted band size: 39kDa. Observed band size: 39kDa.



Immunocytochemistry of NUDT9(green) in HeLa cells using NUDT9 Rabbit pAb at dilution 1/50, and DAPI(blue)