

APOA4 antibody



Catalog Number: 107062

Product name

APOA4 antibody

Immunogen

[Human APOA4 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

APOA4 Mouse Monoclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 45-46 kDa

Preparation

This antibody was obtained by immunization of APOA4 recombinant protein (Accession Number: NM_000482). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107062(APOA4 Antibody) at dilution of 1:4000