

Anti-Bcr Rabbit antibody



Catalog Number: 168062

Product name

Anti-Bcr Rabbit antibody

Specificity

Human

Antibody description

Bcr Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Bcr

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/2000


IHC: 1/20

ICC/IF: 1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of Bcr in
A549,HL-60,U2OS,C6 cell lysates using Bcr Rabbit
mAb(1:1000 diluted).Predicted band
size:143kDa.Observed band size:160kDa.