

# PQBP1 antibody



Catalog Number: 114162

## Product name

PQBP1 antibody

## Immunogen

[Human PQBP1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PQBP1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells, mouse skeletal muscle tissue, rat brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 36-38 kDa

## Preparation

This antibody was obtained by immunization of PQBP1 recombinant protein (Accession Number: XM\_011543884). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114162(PQBP1 antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:114162(PQBP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)