

L2HGDH antibody



Catalog Number: 112262

Product name

L2HGDH antibody

Recommended Dilution:

WB: 1:200-1:2000

Immunogen

IHC: 1:20-1:200

[Human L2HGDH Recombinant protein \(His tag\)](#)


Validations

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

L2HGDH Rabbit Polyclonal antibody. Positive WB detected in mouse small intestine tissue, rat brain tissue, SGC-7901 cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 46-48 kDa




Preparation

This antibody was obtained by immunization of L2HGDH recombinant protein (Accession Number: BC006117). Purification method: Antigen affinity purified.

mouse small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:112262(L2HGDH antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human brain using Catalog No:112262(L2HGDH antibody) at dilution of 1:100 (under 10x lens)

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions



L2HGDH antibody



Catalog Number: 112262

Immunohistochemical of paraffin-embedded

human brain using Catalog No:112262(L2HGDH antibody) at dilution of 1:100 (under 40x lens)