

NUDT6 antibody



Catalog Number: 113362

Product name

NUDT6 antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:50-1:500

Antibody description

NUDT6 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive IF detected in HepG2 cells. Positive IP detected in HepG2 cells. Positive WB detected in Jurkat cells, HepG2 cells, human heart tissue, rat heart tissue, rat liver tissue. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of NUDT6 recombinant protein (Accession Number: NM_007083). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP, IF


Dilutions

Recommended Dilution:

Validations





Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113362(NUDT6 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:113362(NUDT6 antibody) at dilution of 1:50 (under 10x lens)

NUDT6 antibody

Catalog Number: 113362



IP Result of anti-NUDT6 (IP:Catalog No:113362, 4ug; Detection:Catalog No:113362 1:500) with HepG2 cells lysate 4000ug.

Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:113362(NUDT6 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)