

Anti-CACNA1H antibody



Catalog Number: 176362

Product name

Anti-CACNA1H antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to CACNA1H

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cacna1h. aa range:462-511

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1:: Western Blot analysis of various cells using T-type Ca⁺⁺ CP α 1H Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000