

TSPAN33 antibody



Catalog Number: 116362

Product name

TSPAN33 antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TSPAN33 Rabbit Polyclonal antibody. Positive IHC detected in human hepat°Cirrhosis tissue. Positive WB detected in human brain tissue, human blood tissue. Observed molecular weight by Western-blot: 64-72kd

Preparation

This antibody was obtained by immunization of TSPAN33 recombinant protein (Accession Number: NM_178562). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116362(TSPAN33 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:116362(TSPAN33 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog

TSPAN33 antibody



Catalog Number: 116362

No:116362(TSPAN33 antibody) at dilution of 1:50 (under 40x lens)