

Anti-NUP98 (3B8) Mouse antibody



Catalog Number: 168072

Product name

WB, IP

Anti-NUP98 (3B8) Mouse antibody

Dilutions

Specificity

WB: 1/1000

Human

Antibody description

NUP98 (3B8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human NUP98 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Immunoprecipitation analysis of HeLa cell lysates using NUP98 mouse mAb.



Western blot detection of NUP98 in HeLa, Jurkat and MCF7 cell lysates using NUP98 mouse mAb (1:1000 diluted). Predicted band size: 98kDa. Observed band size: 98kDa.