

# FAM133A antibody



Catalog Number: 110472

## Product name

FAM133A antibody

## Immunogen

[Human FAM133A Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

FAM133A Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue.

## Preparation

This antibody was obtained by immunization of FAM133A recombinant protein (Accession Number: XM\_011530926). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, IHC


## Dilutions

Recommended Dilution:  
IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110472(FAM133A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110472(FAM133A Antibody) at dilution of 1:50 (under 40x lens)