

# Anti-AKT Rabbit antibody



Catalog Number: 167472

## Product name

Anti-AKT Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

AKT Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human AKT1/2/3

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions


WB: 1/1000

ICC/IF: 1/200

## Validations



Western blot detection of AKT1/2/3 in LnCap cell lysates using AKT Rabbit pAb(1:500 diluted). Predicted band size: 56kDa. Observed band size: 56kDa.



Immunocytochemistry of AKT1/2/3(green) in Hela cells using AKT Rabbit pAb at dilution 1/50, and DAPI(blue)