

SGOL1 antibody



Catalog Number: 115172

Product name

SGOL1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SGOL1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, COLO 320 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 64kd

Preparation

This antibody was obtained by immunization of SGOL1 recombinant protein (Accession Number: XM_047447487). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115172(SGOL1 antibody) at dilution of 1:400