

# Human FUT10 (His Tag) recombinant protein



Catalog Number: 501972

## General Information

### Gene Name Synonym

Fucosyltransferase X; Galactoside 3-L-fucosyltransferase 10

### Protein Construction

A DNA sequence encoding the luminal domain of human FUT10 (Q6P4F1-1) (Leu 32-Asp 479) was expressed, with a polyhistidine tag at the N-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

His

### Molecular Mass

The recombinant human FUT10 consists of 464

amino acids and has a calculated molecular mass of 54.7 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rh FUT10 is approximately 55-65 kDa due to different glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human FUT10 / Fucosyltransferase 10 Protein (His Tag) SDS-PAGE