

WNT2 Polyclonal Antibody



Catalog Number: 164672

Product name

WNT2 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to WNT2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 26-360 of human WNT2 (NP_003382.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - WNT2 Polyclonal Antibody


Western blot analysis of extracts of SKOV3 cells, using WNT2 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 60s.

Immunohistochemistry - WNT2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using WNT2 antibody at dilution of 1:200



(40x lens).

WNT2 Polyclonal Antibody



Catalog Number: 164672

Immunohistochemistry - WNT2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using WNT2 antibody at dilution of 1:200 (40x lens).