

MYH6 Polyclonal Antibody



Catalog Number: 161372

Product name

MYH6 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to MYH6

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human MYH6 (NP_002462.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:1000 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - MYH6 Polyclonal Antibody

Immunofluorescence analysis of rat heart using MYH6 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.