

# BCRP/ABCG2 antibody



Catalog Number: 107082

## Product name

BCRP/ABCG2 antibody

## Specificity

Human; other species not tested.

## Antibody description

BCRP/ABCG2 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Positive IP detected in NCCIT cell. Observed molecular weight by Western-blot: 70-80; 140-160 kDa

## Preparation

This antibody was obtained by immunization of BCRP/ABCG2 recombinant protein (Accession Number: NM\_004827). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107082(BCRP,ABCG2 antibody) at dilution of



IP Result of anti-BCRP,ABCG2 (IP:Catalog No:107082, 4ug; Detection:Catalog No:107082 1:1000) with NCCIT cell lysate 2000ug.