

PDIA6 antibody



Catalog Number: 113682

Product name

PDIA6 antibody

Recommended Dilution:

WB: 1:1000-1:10000

Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:1000-1:10000

IHC: 1:20-1:200

Antibody description

PDIA6 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, L02 cells, mouse heart tissue, mouse kidney tissue, mouse liver tissue, mouse lung tissue, mouse skin tissue. Positive IP detected in HEK-293 cells. Positive IHC detected in human stomach tissue, human spleen tissue, human stomach cancer tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 54 kDa

IF: 1:20-1:200

Preparation

This antibody was obtained by immunization of PDIA6 recombinant protein (Accession Number: NM_005742). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113682(PDIA6 antibody) at dilution of 1:2000




Immunohistochemical of paraffin-embedded human stomach using Catalog No:113682(PDIA6 antibody) at dilution of 1:50 (under 10x lens)


PDIA6 antibody




Catalog Number: 113682



Immunohistochemical of paraffin-embedded human stomach using Catalog No:113682(PDIA6 antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-PDIA6 (IP:Catalog No:113682, 4ug; Detection:Catalog No:113682 1:3000) with HEK-293 cells lysate 3000ug.



Immunofluorescent analysis of HeLa cells using Catalog No:113682(PDIA6 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)