

CDK13 Polyclonal Antibody



Catalog Number: 160382

Product name

Rabbit IgG

CDK13 Polyclonal Antibody

Applications

Specificity

WB

Human

Dilutions

Antibody description

WB 1:1000 - 1:2000

Polyclonal antibody to CDK13

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 620-900 of human CDK13 (NP_003709.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Validations



Western blot - CDK13 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CDK13 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.