

IL-28A antibody



Catalog Number: 111782

Product name

IL-28A antibody

Specificity

Human; other species not tested.

Antibody description

IL-28A Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of IL-28A recombinant protein (Accession Number: NM_172138). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111782(IL28A Antibody) at dilution of 1:1000