

TULP2 antibody



Catalog Number: 116492

Product name

Polyclonal

TULP2 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

TULP2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, human liver tissue, human lung tissue, Jurkat cells. Observed molecular weight by Western-blot: 59kd

Dilutions


Recommended Dilution:

WB: N/A

Preparation

This antibody was obtained by immunization of TULP2 recombinant protein (Accession Number: NM_003323). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116492(TULP2 antibody) at dilution of 1:300

Clonality