

NOL12 antibody



Catalog Number: 113292

Product name

NOL12 antibody

IP: 1:200-1:1000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NOL12 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HepG2 cells, SH-SY5Y cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of NOL12 recombinant protein (Accession Number: NM_024313). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113292(NOL12 antibody) at dilution of 1:1000



IP Result of anti-NOL12 (IP:Catalog No:113292, 4ug; Detection:Catalog No:113292 1:300) with HeLa cells lysate 4000ug.



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog

NOL12 antibody



Catalog Number: 113292

No:113292(NOL12 Antibody) at dilution of 1:50

and Alexa Fluor 488-conjugated AffiniPure Goat
Anti-Rabbit IgG(H+L)