

# PCDH9 antibody



Catalog Number: 113692

## Product name

PCDH9 antibody

## Immunogen

[Human PCDH9 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PCDH9 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 150-180 kDa

## Preparation

This antibody was obtained by immunization of PCDH9 recombinant protein (Accession Number: NM\_203487). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113692(PCDH9 Antibody) at dilution of 1:600



IP Result of anti-PCDH9 (IP:Catalog No:113692, 4ug; Detection:Catalog No:113692 1:1000) with mouse brain tissue lysate 5000ug.