

USP9X antibody



Catalog Number: 116692

Product name

USP9X antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

USP9X Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, HeLa cells, human brain tissue, K-562 cells. Positive IF detected in HeLa cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 260-290 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001039590). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116692(USP9X antibody) at dilution of 1:500



IP Result of anti-USP9X (IP:Catalog No:116692, 5ug; Detection:Catalog No:116692 1:500) with HEK-293 cells lysate 2400ug.




USP9X antibody




Catalog Number: 116692

Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:116692(USP9X Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded

human pancreas cancer tissue slide using Catalog No:116692(USP9X Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:116692(USP9X Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG