

Anti-TRIM24 Rabbit antibody



Catalog Number: 167692

Product name

Rabbit IgG

Anti-TRIM24 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF

Human, Mouse

Dilutions

Antibody description

WB: 1/1000

TRIM24 Rabbit monoclonal antibody

IHC: 1/50

Preparation

ICC/IF: 1/20

Antigen: Recombinant protein of human TRIM24

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of TRIM24 in HeLa,A549 cell lysates using TRIM24 Rabbit mAb(1:1000 diluted).Predicted band size:117kDa.Observed band size:117kDa.