

C20orf112 antibody



Catalog Number: 108692

Product name

C20orf112 antibody

Immunogen

[Human C20orf112 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

C20orf112 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of C20orf112 recombinant protein (Accession Number: BC065370). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108692(C20orf112 antibody) at dilution of 1:1000