

POLI antibody



Catalog Number: 114003

Product name

POLI antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

POLI Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human heart tissue. Positive IP detected in mouse testis tissue. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 87 kDa

Preparation

This antibody was obtained by immunization of POLI recombinant protein (Accession Number: NM_001351632). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114003(POLI antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:114003(POLI antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded

POLI antibody



Catalog Number: 114003

human pancreas cancer using Catalog No:114003(POLI antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-POLI (IP:Catalog No:114003, 4ug; Detection:Catalog No:114003 1:1000) with mouse testis tissue lysate 4000ug.

