

ZNF549 antibody



Catalog Number: 117003

Product name

ZNF549 antibody

Immunogen

[Human ZNF549 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZNF549 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, A549 cells, HeLa cells, SMMC-7721 cells. Positive IF detected in A549 cells. Observed molecular weight by Western-blot: 74 kDa

Preparation

This antibody was obtained by immunization of ZNF549 recombinant protein (Accession Number: NM_001199295). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:117003(ZNF549 Antibody) at dilution of 1:600



Immunofluorescent analysis of A549 cells using Catalog No:117003(ZNF549 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)