

MRFAP1 antibody



Catalog Number: 112803

Product name

MRFAP1 antibody

Immunogen

[Human MRFAP1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

MRFAP1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 15kd

Preparation

This antibody was obtained by immunization of MRFAP1 recombinant protein (Accession Number: NM_001272053). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112803(MRFAP1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112803(MRFAP1 antibody) at dilution of 1:50 (under 10x lens)