

Anti-PSMC5 Rabbit antibody



Catalog Number: 168803

Product name

Anti-PSMC5 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PSMC5 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human PSMC5

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000-1/10000

ICC/IF: 1/50


FC: 1/20-1/100

IP: 1/20

Validations



Western blot detection of PSMC5 in C6,HeLa cell lysates using PSMC5 Rabbit mAb(1:1000 diluted).Predicted band size:46kDa.Observed band size:46kDa.



Immunocytochemistry of PSMC5(green) in HeLa cells using PSMC5 Rabbit mAb at dilution 1/50, and DAPI(blue)