

RBPJ antibody



Catalog Number: 107503

Product name

RBPJ antibody

Immunogen

[Human RBPJ Recombinant protein \(His tag\)](#)

Specificity

human; other species not tested.

Antibody description

RBPJ Mouse Monoclonal antibody. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of RBPJ recombinant protein (Accession Number: XM_017008174). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107503(RBPJ antibody) at dilution of 1:1000