

Anti-Cytokeratin 18 Rabbit antibody



Catalog Number: 168903

Product name

WB, ICC/IF

Anti-Cytokeratin 18 Rabbit antibody

Dilutions

Specificity

1:2000-1:10000

Human

1:20-1:50

Antibody description

Cytokeratin 18 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Cytokeratin 18

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Rabbit IgG

Applications

Validations



Western blot detection of Cytokeratin 18 in K562, HeLa cell lysates using Cytokeratin 18 Rabbit mAb(1:1000 diluted). Predicted band size: 48kDa. Observed band size: 48kDa.



Immunocytochemistry of Cytokeratin 18(green) in HeLa cells using Cytokeratin 18 Rabbit mAb at dilution 1/200, and DAPI(blue)